# Emissivity Meter TSS-5X-3



### **Typical Applications**

- Easy emissivity measurement at normal room temperature.
- •In material research, designing and planning, thermo-analysys, insulation material research etc.
- ●To ascertain surface condition of material.
- ●in space research, semi-conductor and nuclear energy research.

Technical Specifications	
Measurement system	Reflective infrared energy measurement system
Spectrum range	2 ~ 2 2 μm
Measurement Range	Emissivity: 0.00 ~ 1.00
Accuracy	less than ± 1
Target size	φ 15mm
Measurment distance	12mm (fixed distance from detection head)
Ambient temperature	$10\sim40^\circ\mathrm{C}$ (room temperature)
Measurement Display	LED Digital Display
Output	$0 \sim 0.1 \text{V}, \ 0 \sim 1 \text{V}$ Full scale
Operating temperature/RH	$10 \sim 45^{\circ}\text{C} \cdot 35 \sim 85\%$ (ron-condensiting)
Power supply	AC100 $\sim 240 V_{c}$ 50 $\angle$ 60Hz, 30W (2000L mm cord with connector/plug)
Measuremnt/mass	Detection head: $\phi$ 51 × 137L mm、0.6kg (3000L mm with robot cable)
	Measuring unit 156H × 306W × 230D mm、5.0kg
Accessory	Emissivity standard piece case 84H × 238W × 185D mm、1.1kg
	(Emissivity standard) Emissivity 0.06、0.97 Standard piece One peace each

## Emissivity Meter TSS-5X-3 | Structure


#### Infrared ray radiation source

On heating hemispheric face of full radiator furnace at a constant temperature, it radiates infrared energy to specimen//target and

#### **Infrared detector**

Infrared radiation is reflected by specimen and part of this enter a top hole of hemispheric face of radiator furnace.

Infrared detector detects the energy at constant ratio.



#### **Emissivity Operation Circuit**

Emissivity of the spcimen/target is culculated from outputs of reflected energy detected.

Note: Relational expression of emissivity of specimen/target (γ) and emissivity (  $\epsilon$ ) is :  $\epsilon$ =1 $-\gamma$ 

#### Calibration of displayed emissivity value

Calibration of emissivity of the specimen/target can be done using standard blocks. (  $\varepsilon$  =0.06 and 0.94).

#### **Emissivity Display**

Display emissivity of specimen/target with digital.









DUE TO CONTINUOUS PRODUCT IMPROVEMENT, THE DESIGN AND TECHNICAL SPECIFICATIONS

ARE SUBJECT TO CHANGE WITHOUT PRIOR NOTICE

#### TASHIKA CO., LTD.

1-12, Kaiyo-cho, Ashiya, 659-0035, JAPAN

Tel: +81-797-23-9035 Fax: +81-797-23-2105

e-mail: sales@tashika.co.jp URL: www.tashika.co.jp